

仕 様 書

抗 PreS2 モノクローナル抗体 (マウス)

コード	2APS42
クローン番号	Hyb-5520
由来	HBV の PreS2 領域の合成ペプチドを免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを BALB/c マウスに接種して得られた腹水
特異性	HBV 関連 PreS2 抗原に特異的に反応
精製	硫酸アンモニウム塩析精製
Ig Class	IgG ₁
溶媒	Saline (0.1% NaN ₃ 含有)
蛋白濃度	0.5 mg/mL
容量	0.2 mL
保管	-20 °C以下 (凍結融解の繰り返しはさけてください)
用途	ELISA、Western blotting

参考文献

- 1) Ohnuma H, et al: Gastroenterology **90**: 695-701, 1986.
- 2) Takai E, et al: J Immunol Methods **95**: 23-30, 1986.
- 3) Machida A, et al: Mol Immunol **24**: 523-529, 1987.



〒112-0004 東京都文京区後楽 1-1-10
Tel: 03-3814-4081 Fax: 03-3814-5957
e-mail: info@tokumen.co.jp
<https://www.tokumen.co.jp/>

本製品は研究用에만使用できます。人やその他の動物の診断・治療に使用することはできません。

Specification Sheet

Monoclonal Antibody Against PreS2

Code	2APS42
Clone No.	Hyb-5520
Source	Ascites obtained by inoculating BALB/c mice with hybridoma produced by fusing spleen cells from BALB/c mice immunized with synthetic peptide in the PreS2 region of HBV and mouse myeloma cells (NS-1)
Specificity	This monoclonal antibody specifically reacts with HBV PreS2 antigen.
Purification	Ammonium sulfate precipitation
Ig Class	IgG ₁
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	0.5 mg/mL
Package	0.2 mL
Storage	Below -20°C. Avoid repeated freeze/thaw cycles.
Application	ELISA, Western blotting

- Reference**
- 1) Ohnuma H, et al: Gastroenterology **90**: 695-701, 1986.
 - 2) Takai E, et al: J Immunol Methods **95**: 23-30, 1986.
 - 3) Machida A, et al: Mol Immunol **24**: 523-529, 1987.



INSTITUTE OF IMMUNOLOGY Co., LTD.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004

Tel: +81-3-3814-4081 Fax: +81-3-3814-5957

e-mail: info@tokumen.co.jp

<https://www.tokumen.co.jp/en>

This product is for in vitro research use only, and is not intended for use in the diagnosis or treatment of humans or animals.