Specification Sheet

Purified HBs Antigen

Code 6Z11

Manufacturing

HBsAg positive human plasma

Process

Purification Potassium bromide equilibrium density-gradient centrifugation.

Sucrose density-gradient rate zonal centrifugation.

Gel filtration.

Buffer Saline (containing 0.1% NaN₃)

Protein 1.0 mg/mL

Concentration

Package 1.0 mL

fe

Storage Below -20°C

Notice This product is inactivated.

It is recommended to be handled as if it were potentially infectious,

because it is derived from human plasma.



1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004 Tel: +81-3-3814-4081 Fax: +81-3-3814-5957 e-mail: info@tokumen.co.jp http://www.tokumen.co.jp/english/

6Z111404-01

This product is for *in vitro* research use only, and is not intended to use for diagnosis or treatment for humans nor animals.