Specification Sheet

Purified HBs Antigen

Code	6Z11
Source	HBsAg positive human plasma
Purification	Potassium bromide equilibrium density-gradient centrifugation. Sucrose density-gradient rate zonal centrifugation. Gel filtration.
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	1.0 mg/mL
Package	1.0 mL
Storage	Below -20°C. Avoid repeated freeze/thaw cycles.
Notice	The source sera are tested negative for HCV antibody, HIV-1 antibody, and HIV-2 antibody. Furthermore, this product is inactivated by heat-treatment for 10 hours at 60°C. However it is recommended to be handled as if it were potentially infectious.



1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004 Tel: +81-3-3814-4081 Fax: +81-3-3814-5957 e-mail: info@tokumen.co.jp https://www.tokumen.co.jp/en

This product is for in vitro research use only, and is not intended for use in the diagnosis or treatment of humans or animals.